

DLA - Proficiency Tests GmbH

Hauptstr. 80 23845 Oering/Germany

Tel: +49-(0)4532/9183358 Mob: +49-(0)171/1954375 Fax: +49-(0)4102/9944976 eMail: pt@dla-lvu.de Internet: www.dla-lvu.de

DLA, Hauptstr. 80, 23845 Oering/Germany

Datum / Date: 2. Jul. 2024

Dear participants,

Please find enclosed the material for the proficiency test (PT):

DLA ptSU08 (2024) - Minerals, Trace Elements and Heavy Metals in Turmeric, approx. 20 Elements incl. As, Cd, Hg and Pb

The two portions contain identical samples with the parameters Al, As, Ba, Ca, Cd, Cr, Cu, Fe, Hg, I, K, Mg, Mn, Na, Ni, P, Pb, S, Se and Zn in the matrix of Turmeric (powder). The analysis method is optional.

Please note the attached information on the proficiency test.

<u>New:</u> Please enter your final results online in our <u>PT customer portal</u> **my DLA | participant's portal**. You will receive further information on this by e-mail, in particular about access to the portal.

## Last deadline is August 30th 2024.

After the deadline no results can be accepted.

We are looking forward to any suggestions or questions! We wish you a successful performance of the proficiency test!

Kind regards,

Matthias Besler-Scharf & Alexandra Scharf

On behalf of the DLA-Team

## **DLA - Proficiency Tests GmbH**



Hauptstr. 80 23845 Oering/Germany

Tel: +49-(0)4532/9183358 Mob: +49-(0)171/1954375 Fax: +49-(0)4102/9944976 eMail: pt@dla-lvu.de Internet: www.dla-lvu.de

## Information on the Proficiency Test (PT)

PT number	DLA ptSU08-2024
PT name	Minerals, Trace Elements and Heavy Metals in Turmeric
Sample matrix*	Samples I + II: Turmeric Powder
Number of samples and sample amount	2 identical samples I + II, 8 g each.
Storage	Samples I + II: room temperature
Intentional use	Laboratory use only (quality control samples)
Parameter	quantitative: Al, As, Ba, Ca, Cd, Cr, Cu, Fe, Hg, I, K, Mg, Mn, Na, Ni, P, Pb, S, Se and Zn
Methods of analysis	Analytical methods are optional
Notes to analysis	The analysis of PT samples should be performed like a routine laboratory analysis.  In general we recommend to homogenize a representative sample amount before analysis according to good laboratory practice, especially in case of low sample weights.
Result table	The results for sample I and II as well as the final results calculated as mean of the double determination (samples I and II) should be filled in the result entry table. The recovery rates, if carried out, have to be included in the calculation.
Units	mg/kg
Number of significant digits	at least 2
Further information	For information please specify:  - Date of analysis  - DLA-sample-numbers (for sample I and II)  - Limit of detection  - Assignment incl. Recovery  - Recovery with the same matrix  - Method is accredited
Result submission	online via my DLA   participant's portal (https://my.dla-pt.com) you will receive further information about the access by e-mail
Last Deadline	the latest August 30th 2024
Evaluation report	The evaluation report is expected to be completed 6 weeks after deadline of result submission and sent as PDF file by e-mail.
Coordinator and contact person of PT	Matthias Besler-Scharf PhD

<sup>\*</sup> Control of mixture homogeneity and qualitative testings are carried out by DLA. Any testing of the content, homogeneity and stability of PT parameters is subcontracted by DLA.

DLA - Proficiency Tests GmbH page 2 of 2