

DLA - Proficiency Tests GmbH

Hauptstr. 80 23845 Oering/Germany

Tel: +49-(0)4532/9183358 Mob: +49-(0)171/1954375 Fax: +49-(0)4102/9944976 eMail: pt@dla-lvu.de Internet: www.dla-lvu.de

DLA, Hauptstr. 80, 23845 Oering/Germany

Datum / Date: 9. Sep. 2025

Dear participants,

Please find enclosed the material for the proficiency test (PT):

### ptAU05 (2025) - Steviosides in Instant Tea Drink Powder

There are **two identical samples I and II** with the parameters steviolglycosides (Sum of Steviol Glycosides (as steviol equivalents), Stevioside (as steviol equivalents) and Rebaudioside A (as steviol equivalents)) to be determined. The matrix is instant tea drink powder. All contents should be given as steviol equivalents. The analysis method is optional.

### Please note the attached information on the proficiency test.

Please enter your final results online in our <u>PT customer portal</u> **my DLA | participant's portal**. You will receive further information on this by e-mail, in particular about access to the portal.

#### Last deadline is October 2025-10-24.

After the deadline no results can be accepted.

We are looking forward to any suggestions or questions! We wish you a successful performance of the proficiency test!

Kind regards,

Matthias Besler-Scharf & Alexandra Scharf

On behalf of the DLA-Team

# **DLA - Proficiency Tests GmbH**



Hauptstr. 80 23845 Oering/Germany

Tel: +49-(0)4532/9183358 Mob: +49-(0)171/1954375 Fax: +49-(0)4102/9944976 eMail: pt@dla-lvu.de Internet: www.dla-lvu.de

# Information on the Proficiency Test (PT)

PT number	DLA ptAU05 (2025)
PT name	Steviosides in Instant Tea Drink Powder
Sample matrix*	Samples I + II: Instant tea drink powder (matcha tea with lemon flavor) / Ingredients: Maltodextrin, matcha powder, citric acid, stevia, lemon flavor
Number of samples and sample amount	2 identical samples I + II, 25 g each.
Storage	Samples I + II: room temperature (PT period), cooled 2 - 10°C (long term)
Intentional use	Laboratory use only (quality control samples)
Parameter	quantitative: Stevioside, rebaudioside A and sum of all steviol glycosides (as steviol equivalents)
Methods of analysis	Analytical methods are optional
Notes to analysis	The analysis of PT samples should be performed like a routine laboratory analysis. In general we recommend to homogenize a representative sample amount before analysis according to good laboratory practice, especially in case of low sample weights.
Result table	The results for sample I and II as well as the final results calculated as mean of the double determination (samples I and II) should be filled in the result entry table. The recovery rates, if carried out, have to be included in the calculation.
Units	mg/kg
Number of significant digits	at least 2
Further information	For information please specify:  - Date of analysis  - DLA-sample-numbers (for sample I and II)  - Limit of detection  - Assignment incl. Recovery  - Recovery with the same matrix  - Method is accredited
Result submission	online via my DLA   participant's portal (https://my.dla-pt.com) you will receive further information about the access by e-mail
Last Deadline	the latest October 2025-10-24
Evaluation report	The evaluation report is expected to be completed 6 weeks after deadline of result submission and will be provided as a PDF file in the DLA   Participant Portal (https://my.dla-pt.com/).
Coordinator and contact person of PT	Matthias Besler-Scharf PhD

<sup>\*</sup> Control of mixture homogeneity and qualitative testings are carried out by DLA. Any testing of the content, homogeneity and stability of PT parameters is subcontracted by DLA.

DLA - Proficiency Tests GmbH page 2 of 2